



Multiplex Planar Microarray Assay for Simultaneous Detection of Autoantibodies in ANCA-associated Vasculitis and Anti-GBM or -C1q-associated Glomerulonephritis

Terryberry J¹; Damoiseaux J²; Cervera R³; Kahama A¹; Carmichael S¹; McCabe S¹; Sawyer J;¹ Lea P¹.

1. SQI Diagnostics Toronto, ON Canada 2. Laboratory of Clinical Immunology, University Hospital Maastricht, Netherlands

3. Department of Medicine, University of Barcelona, Spain

Introduction:

Simultaneous measurement of a number of correlated analytes is useful for reasons of reduced specimen volume and reagent costs, and enhanced clinical utility wherein a combination of markers facilitates disease- or organ-oriented profiling, prognostic studies of disease activity and severity and therapeutic monitoring.¹ SQI Diagnostics' IgX PLEXTM planar microarrays allow multiplexing of both antigen and antibody class in an automation compatible 96 well format. Assays which facilitate the detection of multiple biomarkers from multiple patients in a fully automatable platform better reflect the complexity of autoimmune disease processes.² Use of specifically labelled reporter molecules permit simultaneous measurement of autoantibody immunoglobulin classes.³

Objective:

Clinical performance of multiplex testing was evaluated in a novel planar microarray assay (Figure 1) by detecting the frequency of anti-PR3 (proteinase-3), anti-MPO (myeloperoxidase), anti-GBM (glomerular basement membrane antigen alpha3-collagen IV) and anti-C1q IgG antibodies in diagnostic, clinically defined (Chapel Hill and American College of Rheumatology) subsets of vasculitis and nephritis.

Materials and Methods:

Vasculitis sample sets included: 40 ANCA-associated vasculitis (AAV) comprised of Wegener's granulomatosis, microscopic polyangiitis and Churg-Strauss syndrome; 25 systemic lupus erythematosus (SLE) plus vasculitis (SLEV); and 20 other vasculitides (OV) comprised of cutaneous vasculitis, Henoch Schonlein purpura and mixed cryoglobulinemia.

Renal sets included: 23 rapidly progressive (crescentic) glomerulonephritis (RPGN type III); 19 Goodpasture's syndrome (GS, RPGN type I); 71 lupus nephritis (LN) and 20 nephropathy patient samples. 29 healthy controls (HC) were also included. Autoantibodies in patient samples were detected with a research use only version of the IgX PLEXTM Vasculitis and Nephritis assay on a manual adaptation of the fully automated SQiDworksTM Diagnostic Platform.

Clinical sensitivity and specificity were determined with specificity measured against relevant disease controls and healthy controls. Overall correct classification was determined as a measure of differential diagnostic efficiency.

Results:

The clinical performance results support the application of IgX PLEXTM planar microarray technology to differential diagnosis involving ANCA testing in AAV and RPGN type III, GBM antibodies in RPGN type I and C1q antibodies in LN. Sensitivity, specificity, positive and negative predictive values, and overall correct classification for all relevant diagnostic comparisons are listed in Table 1. As a serological test to support clinical symptomatology and histopathology results, the IgX PLEXTM Vasculitis and Nephritis Assay was most efficient for diagnosing AAV, RPGN type III and Goodpasture's syndrome.

Discussion:

The majority of patient sera contained autoantibodies typically associated with their disease state (Table 1. sensitivity values). Multiplex testing for 4 vasculitis- and glomerulonephritis-associated autoantibodies in the IgX PLEXTM microarray detected autoantibody multireactivity that may be overlooked in single antigen assays. The importance of the multiplex assay is best exemplified by the detection of IgG anti-C1q or -GBM antibodies in 55% of RPGN type III but not in AAV. Thus, these kidney antibodies help differentiate the two diseases in the spectrum of ANCA positivity. Additionally, several patients classified as having one disease but exhibiting antibodies typically associated with another disease (e.g. ANCA positive SLEV, LN or GS) may have an overlap syndrome (e.g. ANCA disease⁴) and require additional diagnostic testing. Our results confirm the prevalence of ANCA antibodies in SLE with an increased frequency of ANCA PR3 or MPO antibodies seen in SLE + vasculitis (24%) vs. SLE + nephritis (8.5%).⁵

Substantiating an earlier study⁶, our results demonstrate the coexistence of GBM antibodies and C1q antibodies in a limited number of patients (8.5% of LN), as well as coexistence of GBM antibodies and MPO or PR3 antibodies in 4.7% of AAV and RPGN type III patients. Whether this multireactivity contributes to the etiopathogenesis or severity of the associated disease requires further histopathological and prognostic studies.

Obtaining multiple autoantibody results from one patient sample will enhance clinical decision making based on histopathology results and disease activity scores. Additional biomarkers for inflammation, as well as detection of therapeutic monoclonal antibodies can further expand the clinical utility of the microarray.

Conclusion:

SQI's IgX PLEXTM Vasculitis and Nephritis assay provides a sensitive and specific means of detecting MPO and PR3 antibodies in AAV and RPGN type III, and GBM antibodies in Goodpasture's disease. C1q antibody testing was specific to lupus nephritis. The clinical implementation of this microarray assay will aid physicians in characterizing and monitoring vasculitis and nephritis patients. Further studies are ongoing to determine the clinical utility of IgA antibodies to these antigens.

References:

1. Plebani M, Pittoni M, Celadin M, et al. Autoimmun Rev 2009; 8:238-43
2. Gibson DS, Banha J, Penque D, et al. J Proteomics 2010; 73:1045-60
3. Lea P, Keystone E, Madumba S, et al. Clin Rev Allergy Immunol 2009; Dec 9.
4. Falk RJ and Jennette JC. J Am Soc Nephrol 2010; 21:745-52
5. Dubois' Lupus Erythematosus. DJ Wallace, BH Hahn, eds. Lippincott, Williams & Wilkins, 2007 p334
6. Coremans IE, Daha MR, van der Voort EA, et al. Clin Exp Immunol 1992; 87:256-60

Figure 1. IgX PLEXTM Multiplex Planar Microarray Assay for Vasculitis and Nephritis




Table 1. Clinical Sensitivity, Specificity and Overall Correct Classification of Vasculitis and Nephritis Phenotypes using the IgX PLEXTM Vasculitis and Nephritis Assay

Goodpasture's w/ All other nephritis				AAV w/ SLE+vasculitis and other nephritis			
Anti-GBM	Test Positive	Diagnosis Present	Diagnosis Absent	Anti-MPO or PR3	Test Positive	Diagnosis Present	Diagnosis Absent
	Test Negative	17	6	21		17	7
	Test Negative	2	125	121		8	32
		19	131	150		43	45
	Sensitivity	0.855				Sensitivity	0.805
	Specificity	0.864				Specificity	0.844
	Odds ratio (OR)	177.003				Odds ratio (OR)	21.714
	Likelihood Ratio (+LR)	15.605				Likelihood Ratio (+LR)	5.163
	Likelihood Ratio (-LR)	0.268				Likelihood Ratio (-LR)	0.222
	Positive Predictive Power (PPP)	0.759				Positive Predictive Power (PPP)	0.821
	Negative Predictive Power (NPP)	0.884				Negative Predictive Power (NPP)	0.825
	Overall Correct Classification (OCC)	0.917				Overall Correct Classification (OCC)	0.824
Lupus nephritis w/ All other nephritis				RPGN type III w/ All other nephritis			
Anti-C1q	Test Positive	Diagnosis Present	Diagnosis Absent	Anti-MPO or PR3	Test Positive	Diagnosis Present	Diagnosis Absent
	Test Negative	57	11	68		21	13
	Test Negative	14	48	62		3	37
		71	59	130		24	50
	Sensitivity	0.803				Sensitivity	0.807
	Specificity	0.81				Specificity	0.887
	Odds ratio (OR)	17.065				Odds ratio (OR)	52.231
	Likelihood Ratio (+LR)	4.306				Likelihood Ratio (+LR)	7.404
	Likelihood Ratio (-LR)	0.126				Likelihood Ratio (-LR)	0.107
	Positive Predictive Power (PPP)	0.830				Positive Predictive Power (PPP)	0.818
	Negative Predictive Power (NPP)	0.714				Negative Predictive Power (NPP)	0.970
	Overall Correct Classification (OCC)	0.808				Overall Correct Classification (OCC)	0.881